

For the Attention of the Operating Surgeon:

IMPORTANT INFORMATION

ON THE SYNTHES **ELBOW HINGE FIXATOR, MR CONDITIONAL**

DESCRIPTION

The Synthes Elbow Hinge Fixator consists of two rods, which are interconnected through a riveted joint, allowing a hinge-like movement. The Elbow Hinge Fixator provides stabilization of the elbow via a construct including carbon fiber rods and four Schanz screws, two inserted into the humerus and two inserted into the ulna. The Elbow Hinge Fixator is connected to the carbon fiber rods by utilizing either adjustable or combination clamps. A guide wire is used to determine the anatomical joint axis of the elbow and position the Elbow Hinge Fixator in reference to the axis.

INDICATIONS

The Synthes Elbow Hinge Fixator is intended for supplementary treatment of complex, unstable elbow injuries when early functional stress must be limited due to persistent ligament instability.

The indications for guided joint bridging with external fixators are:

- Delayed treatment of dislocated and stiff elbows
- Chronic, persistent joint instability
- Acute joint instability after complex ligament injuries
- Unstable elbow fractures
- Additional stabilization of post-operative unstable internal fixation

The Elbow Hinge Fixator is compatible with the components of the Synthes Large External Fixation System for adults and with the components of the Synthes Medium External Fixation System for children and small stature adults.

MRI INFORMATION

Synthes Elbow Hinge Fixator devices are labeled MR Conditional according to the terminology specified in ASTM F2503-05, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing of Synthes Elbow Hinge Fixator devices demonstrated that when used in the specific configurations stated in Synthes labeling, are MR Conditional. Representative Synthes Elbow Hinge Fixator devices used in a typical construct included; clamps, rods and various attachments. A patient with a Synthes Elbow Hinge Fixator may be scanned safely after placement of the fixator under the following conditions:

- **Static magnetic field** of 1.5-Tesla when the fixation frame is positioned:
 - 7cm or less from within the outside edge of the bore of the MRI at Normal Operating Mode or;
 - Completely outside of the MRI bore in First Level Controlled Mode
- **Static magnetic field** of 3.0-Tesla when the fixation frame is positioned:
 - 7cm or less from within the outside edge of the bore of the MRI at Normal Operating Mode or;
 - Completely outside of the MRI bore in First Level Controlled Mode
- **Highest spatial gradient magnetic field** of 900-Gauss/cm or less
- **Maximum MR system reported** whole body averaged specific absorption rate (SAR) of 2-W/kg for the Normal Operating Mode and 4 W/kg for the First Level Controlled Mode for 15 minutes of scanning
- **Use only whole body RF transmit coil**, no other transmit coils are allowed, local receive only coils are allowed

Note:

In non-clinical testing, the Synthes external fixation frame was tested in several different configurations. This testing was conducted with the construct position 7cm from within the outside edge of the MRI bore.

- The results showed a maximum observed heating for a wrist fixation frame of 6°C for 1.5T and less than 1°C for 3.0T with a machine reported whole body averaged SAR of 2 W/kg

Patients may be safely scanned in the MRI chamber at the above conditions. Under such conditions, the maximal expected temperature rise is less than 6°C. Because higher in-vivo heating cannot be excluded, close patient monitoring and communication with the patient during the scan is required. Immediately abort the scan if the patient reports burning sensation or pain. To minimize heating, the scan time should be as short as possible, the SAR as low as possible and the device should be as far as possible from the edge of the bore. Temperature rise values obtained were based upon a scan time of 15 minutes.

The above field conditions should be compared with those of the user's MR system in order to determine if the item can safely be brought into the user's MR environment. If placed in the bore of the MR scanner during scanning, Synthes MR Conditional External Fixation devices may have the potential to cause artifact in the diagnostic imaging.

All components of Synthes External Fixation frames must be identified as MR Conditional prior to being placed in or near an MR Environment

Artifact Information:

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of

the Synthes Elbow Hinge Fixator construct and it may be necessary to optimize MR imaging parameters in order to compensate for the presence of the fixation frame.

Representative devices used to assemble a typical Synthes Elbow Hinge Fixator frame have been evaluated in the MRI chamber and worst case artifact information is provided below. Overall, artifacts created by Synthes Elbow Hinge Fixator devices may present issues if the MR imaging area of interest is in or near the area where the fixation frame is located.

- For FFE sequence: Scan duration: 3min, TR 100ms, TE 15ms, flip angle 15° and SE sequence: Scan duration: 4min, TR 500ms, TE 20ms, flip angle 70° radio echo sequence, worst-case artifact will extend approximately 5cm from the device.

WARNING

Do not place any radio frequency (RF) transmit coils over the external fixation frame

STERILIZATION

These devices are offered NONSTERILE only. Resterilization of the device may only be performed if the device has been opened, but not used. Resterilization of the device should not be performed if the device packaging is damaged upon receipt or if the device has been contaminated by bodily fluids.

For detailed cleaning and sterilization instructions, please refer to:

www.depuysynthes.com/hcp/cleaning-sterilization or sterilization instructions, if provided. The parameters are for sterilization of the individual nonsterile device.

In Canada, the cleaning and sterilization instructions will be provided with the Loaner shipments.

CAUTION:

**FEDERAL LAW RESTRICTS THIS DEVICE
TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

Manufactured by:

Synthes USA, LLC
1101 Synthes Avenue
Monument, CO 80132

Synthes GmbH
Luzernstrasse 21
4528 Zuchwil, Switzerland
+41 32 720 40 60

Note: For recognized manufacturer, refer to the product label.