# Medtronic


STRATAMR™ VALVES  
ΚΛΑΠΗ STRATAMR™  
VENTILY STRATAMR™  
STRATAMR™ VENTILER  
STRATAMR™-VENTILE  
ΒΑΛΒΙΔΕΣ STRATAMR™  
VÁLVULAS STRATAMR™  
STRATAMR™-KLAPID  
STRATAMR™-VENTTIILIT  
VALVES STRATAMR™  
VENTILI STRATAMR™  
STRATAMR™ SZELEPEK  
VALVOLE STRATAMR™  
„STRATAMR™“ VOŽTUVAI  
STRATAMR™ VÄRSTI  
STRATAMR™ KLEPPEN  
STRATAMR™-VENTILER  
ZASTAWKI STRATAMR™  
VÁLVULAS STRATAMR™  
VALVE STRATAMR™  
ΚΛΑΠΑΝΗ STRATAMR™  
VENTILY STRATAMR™  
VALVULE STRATAMR™  
STRATAMR™ VENTILI  
STRATAMR™-VENTILER  
STRATAMR™ VALFLERİ

INSTRUCTIONS FOR USE / ИНСТРУКЦИИ ЗА УПОТРЕБА / NÁVOD K POUŽITÍ / BRUGSANVISNING / GEBRAUCHSANWEISUNG / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / INSTRUCCIONES DE USO / KASUTUSJUHEND / KAYTTÖOHJE / MODE D'EMPLOI / UPUTE ZA UPOTREBU / HASZNÁLATI ÚTMUTATÓ / ISTRUZIONI PER L'USO / NAUDOJIMO INSTRUKCIJOS / LIETOŠANAS INSTRUKCIJA / GEBRUIKSAANWIJZING / BRUGSANVISNING / INSTRUKCJA UŻYTKOWANIA / INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCȚIUNI DE UTILIZARE / ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ / NÁVOD NA POUŽITIE / NAVODILO ZA UPORABO / UPUTSTVO ZA UPOTREBU / BRUGSANVISNING / KULLANIM TALİMATLARI

**Rx** only


1

## (1) Small StrataMR Valve


2

## (1) Regular StrataMR Valve

3

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## Description

After surgical implantation, the Medtronic Neurosurgery StrataMR valves provide a non-invasive method to address changing patient needs in the management of hydrocephalus. The valve design enables the physician to adjust the valve pressure/performance level using the Medtronic StrataMR Adjustment Tools.

The StrataMR valves are intended to maintain performance level setting under magnetic resonance imaging (MRI) exposure when scanned in accordance with the MR conditions in the labeling.

The StrataMR valves, small (figure 1) and regular (figure 2), provide the full range of performance levels: 0.5, 1.0, 1.5, 2.0, and 2.5. Performance characteristics of the valve are shown in figures 3, 4, 5a and 5b. Pressure/flow performance data with and without siphon effect is indicated. Each valve is individually tested to ensure conformance with labeled performance characteristics. The currently set performance level is discernible through use of the **Indicator Tool** or by x-ray (for corresponding settings, refer to figure 6).

The dome is designed to allow needle puncture via a 25-gauge or smaller noncoring needle. The StrataMR valve can be flushed in either the distal or proximal direction by percutaneous finger pressure. Occluders are located proximally and distally to the central reservoir on the StrataMR valve regular and only distally on the small model to facilitate selective flushing.

## Indications

The Medtronic StrataMR valves and shunts are designed to provide continuous cerebrospinal fluid (CSF) flow from the ventricles of the brain into the right atrium of the heart or the peritoneal cavity. The design enables the physician to non-invasively adjust valve pressure/performance level pre- and post-implantation by using magnetic adjustment tools without the need for radiographic confirmation.

## Contraindications

Shunting of CSF into the peritoneal cavity or other areas of the body should not be carried out if there is infection in any areas in which the various components of the shunt system will be implanted. This includes infections of the scalp and other skin areas through which the shunt system will traverse, the meninges and cerebral ventricles, the peritoneum and intraperitoneal and retroperitoneal organs, the pleura and the blood stream. CSF shunting is contraindicated if there is infection present in any area of the body. Additionally, shunting into the atrium of patients with congenital heart disease or other serious cardiopulmonary abnormalities is contraindicated.

## Warnings and Precautions

The valve performance level setting should be verified following patient exposure to high magnetic fields. Devices known to contain magnets should be kept away from the immediate valve implant location, as they may have an effect on the performance level setting of the StrataMR valve. All magnets have an exponentially decreasing effect on the valve the further away they are located. Common environmental levels of electromagnetic (radio frequency) radiation generated by security scanners, metal detectors, microwave ovens, mobile telephones, high voltage lines, and transformers should not affect the performance level settings.

An incorrect performance level setting may lead to overdrainage or underdrainage of CSF which may give rise to patient complications.

Radiographic confirmation of the pressure setting is recommended as an alternate method to determine the StrataMR valve performance level setting. Risks associated with radiographic imaging include exposure to low levels of radiation, radiation-induced injuries to the skin and underlying tissues, or the possibility of developing a radiation-induced cancer. Care should be taken to use the least amount of radiation exposure needed to produce the image.

When securing catheters to connectors, the encircling ligatures should be securely, but not too tightly fastened, to avoid cutting through the silicone tubing over time.

While manipulating the ventricular catheter at the valve inlet, ensure the distal end does not become lost within the brain tissue to avoid complications.

Valve function and performance level setting should be checked in the event that the valve is subjected to significant mechanical shock or trauma.

Implanting the valve under more than 1 cm of tissue can interfere with subsequent reading and setting of the valve performance level.

The appropriate product and size must be chosen for the specific patient's needs, based on diagnostic tests and physician experience.

Failure to do so may lead to patient complications.

Use Medtronic Neurosurgery StrataMR Adjustment Tools to change the performance level in the StrataMR valve. Refer to StrataMR Adjustment Tools insert for instructions, warnings, precautions and complications. Use of other adjustment tools may result in inaccurate readings, an inability to adjust the StrataMR valve, or may damage the valve.

The performance characteristics of this device may be altered if components or devices of other manufacturers are used in conjunction with this device.

Overdrainage of CSF may predispose development of a subdural hematoma or hydroma, or collapse of the lateral ventricular walls leading to obstruction of the ventricular catheter.

Shunt obstruction may occur in any of the components of the shunt system. The system may become occluded internally due to tissue fragments, blood clots, tumor cell aggregates, bacterial colonization, or other debris. Catheters which contact internal body structures can become kinked or blocked at their tips. Occluded catheters can lead to underdrainage.

Shunt disconnection may also occur due to growth of an infant or child, or physical activities which result in disconnection of the shunt components or withdrawal of a distal catheter from its intended drainage site. Disconnection may cause complications. Disconnected shunt components may migrate.

Malfunction or obstruction of the shunt system may lead to signs and symptoms of increased intracranial pressure if the hydrocephalus is not compensated. In the infant, the common symptoms are increased tension of the anterior fontanelle, congestion of scalp veins, listlessness, drowsiness and irritability, vomiting and nuchal rigidity. In older children and adults, the common symptoms are headaches, vomiting, blurring of vision, nuchal rigidity, deterioration of consciousness, and various abnormal neurological findings.

Patients with CSF shunt systems must be kept under close observation in the postoperative period for signs and symptoms that suggest shunt malfunction. The clinical findings may indicate shunt infection or other complications such as CSF leaks due to the puncture of the dura, shunt obstruction, underdrainage or overdrainage of CSF.

Shunt systems may fail due to mechanical malfunction leading to underdrainage or overdrainage.

Improper use of instruments in the handling or implantation of shunt products may result in the cutting, slitting or crushing of components. Such damage may lead to loss of shunt integrity and necessitate premature surgical revision of the shunt system.

Care must be taken to ensure that particulate or other contaminants (such as latex) are not introduced into shunt components during preimplantation testing or handling. Introduction of contaminants could result in improper performance of the shunt system, foreign body reactions, or allergic reactions. Particulate matter that enters the shunt system may result in shunt occlusion, or may also hold pressure/flow controlling mechanisms open, resulting in overdrainage.

## MRI Safety Information

Non-clinical testing has demonstrated that the StrataMR valves and shunts are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 or 3.0 tesla
- Maximum spatial gradient magnetic field of 1900 gauss/cm (19 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- The valve setting should be verified following MRI exposure

Under the scan conditions defined above, the StrataMR valves and shunts are expected to produce a maximum temperature rise of less than 3° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 30 mm from the StrataMR valve when imaged with a gradient echo pulse sequence and a 3.0 tesla MRI system.

## Optional Preimplantation Performance Level Test

Medtronic does not recommend functional testing of the valve prior to implantation, due to additional handling and potential contamination. If testing is considered necessary, instructions are available upon request.

## Instructions for Use

### Setting the Performance Level before Implantation

The physician must determine the initial performance level setting for each patient and adjust the valve accordingly prior to implantation. StrataMR valve PL settings are adjusted with the StrataMR Adjustment Tools. Refer to the Instructions for Use for the StrataMR Adjustment Tools.

**Warning: StrataMR valves may not be resterilized. It is essential to set the performance level before implantation with the valve sealed in its sterile package.**

### Patency Check

If pre-fill or patency check is desired, place the inlet connector of the valve into filtered, sterile isotonic saline. Depress and release the valve dome repeatedly until fluid flows from the outlet connector (figure 7). To avoid retrograde flow, occlude the valve inlet when depressing the dome. If fluid flows from the outlet connector each time the dome is depressed, the valve is patent.

Testing shows that applying a water column of up to 100 cm positive pressure or 50 cm negative pressure to the shunt system does not impair performance.

**Caution: Do not attempt to aspirate fluid through the StrataMR valve from the distal (outlet) end. Damage to the valve may result.**

**Caution: Particulate matter in solutions used to test valves may result in improper product performance. Take care to maintain sterility and to avoid particulate contamination.**

### Surgical Technique

A variety of surgical techniques may be used in placing the StrataMR valve. The valve is implanted with the flat surface adjacent to the pericranium. Site of placement is at the discretion of the surgeon. To reduce the possibility of post-operative valve migration (e.g., as the result of magnetic influence due to MRI), suture the valve to adjacent tissue by passing a suture through the polyester-fabric-reinforced flanges.

It is suggested that the StrataMR valve be placed in a surgically created loose subgaleal pocket, avoiding compression by the overlying scalp, and not under the scalp incision. The StrataMR valve should not be placed under the skin of the neck, chest, or abdomen.

Placement of a StrataMR valve will have an impact on overall shunt performance (figure 8). The foramen of Monro and the tip of the ventricular catheter are established zero level reference points for intracranial pressure monitoring. Placement of the valve above these reference points will result in an overall increase in the resistance to flow of the shunt system, with potential underdrainage. Placement of the valve below these reference points will result in an overall decrease in the resistance to flow of the shunt system, with potential overdrainage. These factors must be taken into account for proper patient therapy.

**Caution: Valve placement can have significant effects on the performance of the StrataMR valve.**

The resistance properties of StrataMR valves as indicated in the Performance Characteristics charts in this product labeling are specified at the zero reference point.

Connect valve to catheters by inserting integral connectors into catheters. The connectors should be completely covered by catheter tubing. Secure catheters to connectors with encircling ligatures.

Prior to closure, patency of the shunt system should be verified, through observation of flow.

### Puncturing the Valve Dome

The StrataMR valve is designed to allow CSF sampling through the dome by use of a 25-gauge or smaller noncoring needle (figure 9).

The needle should be inserted at an angle of approximately 30° to 45° from the scalp or base of the valve. If the valve will be punctured several times, it is recommended that the needle be inserted at various locations to avoid multiple punctures at a single point. The catheter tubing, occluders, and siphon control diaphragms should **not** be used as injection sites.

**Caution: Low tear strength is a characteristic of most silicone elastomer materials. Care must be taken on insertion and removal of the needle.**

**Caution: The siphon control portion of this device is not designed to allow needle penetration. Puncture of the diaphragms may affect the performance characteristics of the valve and compromise the shunt system.**

### Valve Flushing

The StrataMR valve can be flushed in both the proximal and distal directions by percutaneous depression of the valve dome (figure 10). To selectively flush the StrataMR valve regular, depress and occlude either proximal or distal occluder sections of the valve by percutaneous finger pressure and then depress the valve dome. For the small valve, distal flushing can be achieved by depressing the ventricular catheter tubing for occlusion.

The selective flushing will cause fluid to flow in the direction opposite the occluded side of the valve. If there is noticeable resistance to compression, the catheter being flushed may be occluded.

**Caution: Shunt obstruction may occur in any component of a shunt system and should be diagnosed by clinical findings and diagnostic testing. Valve flushing characteristics may not be adequate to diagnose occlusion of catheters.**

### Setting the Performance Level After Implantation

StrataMR valve pressure/performance level settings are adjusted with the StrataMR Adjustment Tools. Refer to the Instructions for Use which accompany the StrataMR Adjustment Tools for information regarding performance level adjustments after implantation. More than one attempt may be required to read or adjust the valve setting.

**Caution: Excessive swelling or bandages may make it difficult to determine a valve setting. Wait until swelling is reduced, or confirm with x-ray. In addition, scalp thickness can make it difficult to determine valve setting. If this is the case, confirm valve setting with x-ray.**

## Patient Card

Each StrataMR valve is packaged with a patient card, which includes important MR imaging information and a set of adhesive-backed labels. One label for each shunt component can be used for the chart in the operating room, the patient's clinic chart, and the patient card. **The card should be given to the patient to be carried at all times to provide a current patient record of implanted devices.** Record the performance level setting on the patient card during the initial setting and for subsequent revisions to the performance level setting.

## How Supplied

Medtronic Neurosurgery StrataMR valves are packaged **STERILE** and **NON-PYROGENIC** and are intended for **single** (one-time) **use only**.

The enclosed product is designed for single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Do not use if package has been previously opened or damaged. Medtronic Neurosurgery is not responsible for the performance of any product that has been resterilized.

## Patient Education

It is the responsibility of the physician to educate the patient and/or their representative(s) regarding CSF shunting. This should include a description of the complications associated with implantable shunt systems, and an explanation of potential alternative products and treatments.

## Complications

Complications associated with ventriculoperitoneal CSF shunting procedures may be similar to those experienced in any surgical procedure carried out under local and/or general anesthesia including reactions to drugs and anesthetic agents, electrolyte imbalance, and excessive blood loss, particularly in infants.

In CSF shunting procedures, the most common complication is obstruction of the shunt system. Obstruction by brain fragments, blood clots, and/or tumor cell aggregates may occur in any component of the system and lead to underdrainage of the CSF. Obstruction may also occur because of separation of the system components or kinking/coiling of the catheter. This may predispose migration of the ventricular catheter into the lateral ventricle and the distal catheter into the peritoneum or other structure in which the catheter is implanted. Obstruction due to CSF overdrainage may also occur; excessive reduction of CSF pressure can cause the ventricular walls to impinge on the inlet holes in the catheter.

Overdrainage can also predispose the development of a subdural hematoma or hygroma. In the infant, an excessive pressure reduction from overdrainage can cause marked depression of the anterior fontanelle and overriding cranial bones and may convert communicating into obstructive hydrocephalus.

There are other potentially serious complications. Local and systemic infections due to organisms inhabiting the skin, particularly *Staphylococcus epidermidis*, may occur. Other pathogens circulating in the bloodstream may colonize the shunt and, in the majority of patients, require its removal. Rigorous perioperative control of the operating environment and the use of antibiotics at the physician's discretion may avert infection occurrence.

Growth of the infant or child may cause the distal catheter to be withdrawn from the atrium into the internal jugular vein or from the peritoneum into tissue planes where the fluid cannot be absorbed.

## Returned Goods Policy

Products must be returned in unopened packages, with manufacturer's seals intact, to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling. Determination of a product defect or mislabeling will be made by Medtronic Neurosurgery, which determination will be final. Products will not be accepted for replacement or credit if they have been in the possession of the customer for more than 90 days.

## Warranty

**A. Standard Limited Warranty.** Medtronic Neurosurgery warrants to the original end user purchaser ("Purchaser") that Medtronic Neurosurgery's product ("Product") purchased by Purchaser, at the time of delivery to Purchaser, shall be substantially free from defects in material and workmanship. Medtronic Neurosurgery makes no warranty (express, implied, or statutory) for Products that are modified (except as expressly contemplated herein) or subjected to unusual physical stress, misuse, improper operation, neglect, improper testing, use in combination with other products or components other than those for which the Products were designed, or use in any manner or medical procedure for which the Products are not indicated.

**B. Remedy.** Purchaser's exclusive remedy and Medtronic Neurosurgery's sole liability for breach of the foregoing warranty shall be, at Medtronic Neurosurgery's sole option and election, to replace the Product or credit Purchaser for the net amount actually paid for any such Product; provided that (i) Medtronic Neurosurgery is notified in writing within ninety (90) days after Purchaser's receipt of the Product that such Product failed to conform, including a detailed explanation in English of any alleged nonconformity; (ii) such Product is returned to Medtronic Neurosurgery within ninety (90) days after Purchaser's receipt of the Product F.O.B. 125 Cremona Drive, Goleta, California 93117, U.S.A. or as otherwise designated by Medtronic Neurosurgery; and (iii) Medtronic Neurosurgery is reasonably satisfied that the claimed nonconformities actually exist. Except as expressly provided in this paragraph, Purchaser shall not have the right to return Products to Medtronic Neurosurgery without Medtronic Neurosurgery's prior written consent.

**C. Exclusion of Other Warranties. EXCEPT FOR THE LIMITED WARRANTY PROVIDED IN (A) ABOVE, MEDTRONIC NEUROSURGERY GRANTS NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, AND MANUFACTURER SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC NEUROSURGERY NEITHER ASSUMES NOR AUTHORIZES ANY OTHER PERSON TO ASSUME ANY OTHER LIABILITIES ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF ANY PRODUCT.**

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© 2016 Medtronic, Inc.
Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic, Inc.