Instructions for Use

154-7581-1

PALMAZ® GENESIS™ Peripheral Stent

STERILE. The PALMAZ GENESIS Peripheral Stent is sterilized with gamma radiation. Nonpyrogenic. The stent is radiopaque. For single use only. Do not autoclave.

NON-STERILE. The Cordis introducer tube, provided separately, is supplied NON-STERILE and must be sterilized prior to use. The introducer tube should be sterilized by autoclaving in accordance with hospital procedures.

Caution: Not for sale in the USA.

I. Device Name

The device brand name is the PALMAZ GENESIS Peripheral Stent.

II. Description

The PALMAZ GENESIS Peripheral Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is supplied in two nominal lengths: 29 mm and 39 mm. Specific expansion ranges for each stent length are shown in Table 1. The PALMAZ GENESIS Peripheral Stent is sold unmounted for use with all catheter lengths of the Cordis OPTA™ PRO Percutaneous Transluminal Angioplasty (PTA) Balloon Catheter [a .035" (.89 mm) diameter guidewire compatible system]. For further information, see Table 1.

	Table 1											
	ent ription	Ste Leng		Recommended Co Balloon Cath		Infla	ninal ation ssure	Recommend Catheter Introducer ¹ Guiding C	Sheath (CSI) and			
Product Code	Nominal Dia x L (mm)	Unexpanded (mm)	Expanded ³ (mm)	Catheter	Catalog Number ⁴	atm	kPa	CSI Size French (mm)	Guiding Catheter Size French (mm)			
PG295P	5 X 29 6 X 29 7 X 29 8 X 29 9 X 29	29	28 28 27 26 25	OPTA PRO	419-5030* 419-6030* 419-7030* 419-8030* 419-9030*	6 6 6 6	608 608 608 608 608	7F (2.3) 7F (2.3) 7F (2.3) 7F (2.3) 7.5F (2.5)	8F (2.7) 9F (3.0) 9F (3.0) 9F (3.0) 10F (3.3)			
PG395P	5 X 39 6 X 39 7 X 39 8 X 39 9 X 39	39	38 38 37 36 35	OPTA PRO	419-5040* 419-6040* 419-7040* 419-8040* 419-9040*	6 6 6 6	608 608 608 608	7F (2.3) 7F (2.3) 7F (2.3) 7F (2.3) 7.5F (2.5)	8F (2.7) 9F (3.0) 9F (3.0) 9F (3.0) 10F (3.3)			

- 1 Use Cordis introducer tube INTR6 for all product codes and lengths.
- 2 Use of guiding catheters other than those recommended must have an inner lumen diameter of .088" (2.2 mm) or larger for 8F (2.7 mm) guides, .098" (2.5 mm) or larger for 9F (3.0 mm) guides, and .110" (2.8 mm) or larger for 10F (3.3 mm) guides.
- 3 These data are based on *in-vitro* testing.
- 4 Catheter suffix (*) refers to the usable catheter length. **OPTA PRO** is offered in the following lengths: T (40 cm), V (65 cm), S (80 cm), L (110 cm), X (135 cm).

Note: Not all balloon sizes are sold in all catheter lengths.

III. Indications for Use

The **PALMAZ GENESIS** Peripheral Stent is indicated for use in the treatment of atherosclerotic disease of peripheral arteries below the aortic arch and for palliation of malignant neoplasms in the biliary tree.

IV. Contraindications

Peripheral Artery Stent Implantation

Generally, contraindications to percutaneous transluminal angioplasty (PTA) are also contraindications for stent placement. Contraindications include, but may not be limited to:

- · Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- · A lesion that is within an aneurysm or an aneurysm with a proximal or distal segment adjacent to the lesion.

Transhepatic Biliary Stent Implantation

Contraindications associated with the use of transhepatic biliary endoprostheses include:

- · Stenting of a perforated duct where leakage from the duct could be exacerbated by the prosthesis.
- Patients with uncorrected bleeding disorders.
- · Severe ascities.

V. Warnings

- Do not use if the inner package is opened or damaged.
- The PALMAZ GENESIS Peripheral Stent is intended for single use only. Do not resterilize and/or reuse the device.
- Store in a cool, dark, dry place.
- Use the stent prior to the "Use By" date specified on the package.
- Persons with allergic reactions to stainless steel or its components (for example, nickel) may suffer an allergic response to this implant.
- Once fully deployed, the stent cannot be repositioned.
- Stenting across a bifurcation could compromise future diagnostic or therapeutic procedures.
- Avoid stent placement which would obstruct access to a vital side branch.
- It is recommended that stents not be implanted in patients with confirmed pregnancies.

For Peripheral Artery Stent Implantation Only:

- The safety and efficacy of the PALMAZ GENESIS Peripheral Stent have not been established for use in patients for the treatment of atherosclerotic disease of the femoral artery and other superficial arteries.
- As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or the retroperitoneum.
- The stent may cause a thrombus, distal embolization or may migrate from the site of the implant down the arterial lumen.
 Overstretching of the artery may result in rupture and life threatening bleeding.
- · When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.

VI. Precautions

- The device should only be used by physicians who are trained in interventional techniques such as percutaneous transluminal angioplasty, placement of stents, and transhepatic access.
- · The delivery system is not designed for use with power injection systems.
- · Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.
- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception, stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.
- · To assure full expansion, inflate to at least the nominal pressure recommended on the catheter label.
- · Do not exceed the rated burst pressure recommended on the catheter label.
- Ensure stent is securely crimped onto balloon prior to use.
 Do not attempt to remove or readjust the stent once crimped onto the delivery system.
- · When catheters are in the body, they should be manipulated only under fluoroscopy.

For Peripheral Artery Stent Implantation Only:

- When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the
 proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and
 reduces the chance for dislogging the proximal stent.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infections, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.
- Caution should be taken when stenting patients with poor renal function who, in the physician's opinion, may be at risk for a
 contrast medium reaction.
- In patients requiring the use of antacids and/or H2-antagonists before or immediately after stent placement, oral absorption
 of antiplatelet agents (e.g. aspirin) may be adversely affected.

VII. Potential Complications

Potential complications associated with **peripheral artery and transhepatic biliary stent** implantation may include, but are not limited to, the following:

- Sepsis/infection
- · Stent migration/embolization
- · Drug reaction, allergic reaction to contrast medium

Potential complications associated with peripheral artery stent implantation may include, but are not limited to, the following:

- · Embolization of atherosclerotic or thrombotic material
- Acute or sub-acute stent thrombosis
- Amputation
- Arteriovenous fistula
- · Emergency surgery to correct vascular complications
- GI bleeding from anticoagulation/antiplatelet medication
- Hemorrhage/hematoma
- · Artery injury, including perforation and dissection
- Rupture of retroperitoneum or neighboring organ
- Intimal tear
- Pseudoaneurysm formation
- Restenosis of the stented artery
- Thrombus
- Tissue necrosis
- Total occlusion

Potential complications associated with **transhepatic biliary endoprostheses** implantation may include, but are not limited to, the following:

- Stent obstruction secondary to tumor growth through the stent
- Tumor overgrowth at the stent ends
- Bile duct occlusion/obstruction
- Pancreatitis
- Bile duct perforation
- Parenchymal hemorrhage
- Hemobilia
 Peritonitis
- * Ferriorius
- Abscess
- · Rupture, overstretching of duct
- Cholangitis

VIII. Directions for Use

Procedure

Pre-Procedure for Peripheral Artery Stenting

The patient may be started on enteric coated or nonenteric-coated aspirin 81-325 mg daily, one or two days prior to the procedure if deemed appropriate by the physician.

The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

Procedure for Peripheral Artery and Transhepatic Biliary Stenting

- 1. Perform standard diagnostic angiography or percutaneous cholangiogram, as applicable, to evaluate the lesion.
- Select Stent Size
 - a. Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent (see Table 1).

Note: Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.

- Measure the diameter of the reference vessel to determine the appropriate size stent and delivery system.
- 3. Preparation of Stent
 - a. Open the carton to reveal the pouch containing the stent.
 - b. Inspect the stent package for damage to the sterile barrier. Remove the stent from the package and rinse in sterile heparinized saline.
- 4. Preparation of Stent Delivery Catheter See Instructions for Use provided with the recommended balloon catheter.
 - a. Remove the balloon protector.
 - **b.** Attach a stopcock to the catheter's inflation port.
 - c. Open the stopcock and induce negative pressure.

d. Inflate the balloon by holding the proximal end of the catheter above the distal end. Hold the balloon vertically, with the balloon tip pointing down.

Warning: Inflation at a high rate may damage the balloon.

- e. Continue pointing the distal tip of the balloon downward, and deflate the balloon again by inducing negative pressure.
- f. Ensure that the air in the balloon and the inflation lumen is removed. Repeat steps d-f.
- g. While holding negative pressure, close the stopcock. Maintain negative pressure until the stent is mounted.
- h. Looking proximal to distal, manually refold the balloon clockwise around the catheter.
- 5. Preparation of Stent Delivery System
 - a. Slide the stent over the distal end of the balloon, maintaining the balloon fold, until the radiopaque marker bands are equidistant from the ends of the stent.
 - b. Crimp the stent onto the balloon using the thumb and index finger of both hands. Apply even, concentric pressure along the length of the stent while fully rotating the assembly. Ensure that the stent is securely adhered to the balloon.
 - c. Open the stopcock allowing the inflation lumen and balloon to fill with diluted contrast medium.
 - d. Induce negative pressure to check the balloon for leaks, and then slowly release the negative pressure.
 - e. Remove the stopcock from the catheter balloon inflation port.
- 6. Stent Deployment Procedure for Peripheral Artery and Transhepatic Biliary Stenting

Insertion of Cordis Catheter Sheath Introducer (CSI) and Guidewire

- a. Gain access at the appropriate site utilizing the recommended CSI size (see Table 1).
- Caution: Always use a CSI for the implant procedure to protect the puncture site and, in the case of biliary stenting, the liver tract, and to avoid dislodging the stent from the balloon.
- b. Insert the appropriately sized guidewire that is at least twice the length of the delivery system across the lesion to be stented through the CSI.
- **Note:** If the physician determines that predilation is necessary, PTA techniques may be used. Maintain lesion access with the quidewire.

Caution: Avoid overexpansion of the balloon to prevent bleeding, dissection or patient discomfort.

c. Advance the CSI completely across the lesion.

Optional: Select and insert the appropriate size guiding catheter if the CSI utilized during the procedure cannot be advanced to the lesion site

Introduction of Stent Delivery System

- a. Flush guidewire lumen of the delivery system. Position the flared end of the introducer tube over the distal tip of the stented balloon and slide forward to protect the stent during CSI insertion. Backload onto the guidewire. Place the assembly through the CSI hemostatic valve until resistance is met. Carefully advance the stent and delivery system through the introducer tube and hemostatic valve. When the stent has passed into the body of the CSI, remove the introducer tube.
- b. Continue to advance the device over the guidewire through the hemostatic valve and CSI.

Stent Denloymen

- Keeping the CSI immobile, observe fluoroscopically as the stent is advanced through the CSI to the site of the lesion.
 Caution: If excessive resistance is encountered, do not force passage.
- b. Under fluoroscopy, use the balloon marker bands and the radiopaque stent to position the stent centrally within the lesion. During positioning, verify that the stent is still centered within the balloon marker bands and has not been disloded.
- c. Prior to stent expansion, ensure that stent and balloon are completely exposed from the CSI or guiding catheter. Caution: Never re-advance a CSI or guiding catheter over an exposed stent to avoid dislodging the stent.
- d. Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the catheter label. Expand the diameter of the stent to the diameter of the reference vessel.

Caution: Do not exceed the rated burst pressure shown on the catheter label.

Note: It is strongly recommended that the guidewire, the delivery system, or both, remain across the lesion until the procedure is completed and the delivery system is removed from the vessel.

Note: Never use air or any gaseous medium to inflate the balloon.

Note: For peripheral artery stenting, if the patient has significant femoral artery disease or if additional surgical intervention is contemplated, the stent may be placed in an operating room following surgical cut-down and isolation of the femoral artery.

7. Delivery System Withdrawal

- a. After deploying the stent, deflate the balloon by pulling a vacuum, allowing adequate time for the balloon to fully deflate prior to removal
- b. Carefully rotate balloon counterclockwise to ensure separation of the balloon from the stent.
- c. While maintaining negative pressure on the balloon, slowly withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.
- d. Withdraw the deflated delivery system into the CSI or guiding catheter.
- e. Remove and discard the delivery system, guidewire and CSI or guiding catheter from the body.

IX. Magnetic Resonance Imaging (MRI) Compatibility

The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the stainless steel stent should not be greater than those caused by metal surgical clips. In a study of larger size stainless steel stents, no stent migration was observed in an MRI field strength up to 1.5 Tesla.⁵ However, an MRI scan should not be performed until the stent implantation site has had a chance to heal (estimated to be 8 weeks), in order to further minimize the risk of migration.

Protected under one or more of the following U.S. patents: 4,733,665, 4,739,762, 5,102,417, 5,643,312, 5,697,971, 5,879,370, 5,902,332, and other U.S. and/or foreign patents pending.

PALMAZ, PALMAZ GENESIS, and OPTA are trademarks of Cordis Corporation.

5Data on file.

X. DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE CORDIS PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL CORDIS BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CORDIS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Symbols used on product labels and box:

REF Catalog No.

LOT

Lot No.

Use By

For one use only

Sterilized by irradiation

Attention, see Instructions for Use

Cordis Sales / Marketing Offices:

Austria:

Johnson & Johnson Medical Products GmbH Cordis Division Gunoldstrasse 16 A-1190 Wien Telephone 01-360 25-0

Belgium:

Cordis, a Johnson & Johnson Company Johnson & Johnson Medical N.V./S.A. Eikelenbergstraat 20 B-1700 Dilbeek Telephone 02-481 74 00

Canada:

Johnson & Johnson Medical Products 200 Whitehall Drive Markham, Ontario Canada L3R 0T5 Telephone 905-946-1611

European HQ:

Cordis, a Johnson & Johnson Company Johnson & Johnson Medical N.V./S.A. Waterloo Office Park, Building H Drève Richelle 161 B-1410 Waterloo Belgium Telephone 02-352 14 11

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Hong Kong:

Johnson & Johnson Hong Kong, Ltd. Medical Division Room 1816-1819, 18/F Grand Century Place, Tower 1 193, Prince Edward Road West Mongkok, Kowloon Telephone 2738 2818

Italy:

Cordis Italia S.p.A. Via Chiese, 74 I-20126 Milano Telephone 02-64410.1

Japan:

Cordis Japan Johnson & Johnson K.K. East 21 Tower 10th Floor 6-3-2 Toyo, Koto-ku Tokyo 135-0016 Telephone 03-5632-7200

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Canada:

Johnson & Johnson Medical Products 200 Whitehall Drive Markham Ontario Canada L3R 0T5 Telephone 905-946-1611

European HQ:

Cordis, a Johnson & Johnson Company, Johnson & Johnson Medical N.V./S.A., Waterloo Office Park, Building H Drève Richelle 161 B-1410 Waterloo Belgium Telephone 02-352 14 11

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United Kingdom:

Cordis, a Johnson & Johnson Company Johnson & Johnson Medical Ltd. Coronation Road, South Ascot Berkshire SL5 9EY Telephone 01344 871000

USA:

Cordis Corporation P.O. Box 025700 Miami, FL 33102-5700 Telephone 305-824-2000

Cordis Corporation P.O. Box 4917 Warren, NJ 07059-0917 Telephone 908-755-8300

Cordis Operations:

The Netherlands:

Cordis Europa N.V. Oosteinde 8 NL-9301 LJ Roden Telephone 050-5022222

USA:

Cordis Corporation P.O. Box 025700 Miami, FL 33102-5700 Telephone 305-824-2000

EU Authorized Representative Cordis Europa, N.V.

Oosteinde 8 NL-9301 LJ Roden The Netherlands

Manufactured in: USA

RECYCLED 100% Recycled Fibers Including 20% Post Consumer Waste

年9月 2001 © Cordis Corporation 2001 154-7581-1

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Johnson & Johnson Medical BV Postbus 188 NL-3800 AB Amersfoort Telephone 033-450 0729
Johnson & Johnson S.A. Paseo de las doce Estrellas, 5-7 Campo de las Naciones E-28042 Madrid Telephone 91 722 8000
Johnson & Johnson AB Staffans väg 2 SE-191 84 Sollentuna Telephone 08-626 22 00
Johnson & Johnson AG Cordis Division Rotzenbühlstrasse 55 CH-8957 Spreitenbach Telephone 056-417 3207