Medtronic

STRATA II™ VALVE

STRATA II™ VENTIL

STRATA II™-KLEP

VALVE STRATA II™

STRATA II™ VENTIL

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

Description

After surgical implantation, the Medtronic Neurosurgery Strata II valve provides a non-invasive method to address changing patient needs in the management of hydrocephalus. The unique valve design allows the physician to adjust valve pressure/performance level by using the Medtronic Neurosurgery Strata II Adjustment Kit or the StrataVarius Adjustment System. The design helps ensure that the valve's performance level resists inadvertent changes.

The Strata II valve incorporates a ball and cone pressure valve in series with a normally closed siphon control mechanism (Delta® chamber). This combination enables the valve to maintain intraventricular pressure (IVP) within a normal physiological range, regardless of a patient's cerebrospinal fluid (CSF) flow requirements or body position. Flow control is accomplished by combined resistance of the ball and cone and siphon control diaphragm. The degree of resistance determines the performance characteristics of the valve. Retrograde flow is prevented by the ball and cone. The siphon effect is controlled by the two normally closed silicone elastomer diaphragms of the Delta chamber. The diaphragms open in response to positive upstream IVP. The Strata II valves, small and regular, provide the full range of Performance Levels: 0.5, 1.0, 1.5, 2.0, and 2.5. Performance characteristics of the valve are shown in Figures 1 through 2b. Pressure/flow performance data with and without siphon effect is indicated. Each valve is individually tested to ensure conformance with labeled performance characteristics.


1 Strata II Valve Performance at 0 cm H₂O Distal Hydrostatic Pressure

Every Medtronic Neurosurgery valve is tested for its performance characteristics using a water-filled system at the time of manufacture. All valves must conform to labeling specifications in these tests. Subsequent testing of the valve may give different results depending on test conditions and condition of the valve. Testing at the time of manufacture serves to assure consistent, high quality of every Medtronic Neurosurgery valve.


2 Strata II Valve Performance at -50 cm H₂O Distal Hydrostatic Pressure

Every Medtronic Neurosurgery valve is tested for its performance characteristics using a water-filled system at the time of manufacture. All valves must conform to labeling specifications in these tests. Subsequent testing of the valve may give different results depending on test conditions and condition of the valve. Testing at the time of manufacture serves to assure consistent, high quality of every Medtronic Neurosurgery valve.

2a Strata II Valve Performance at 0 cm H₂O Distal Hydrostatic Pressure

2b Strata II Valve Performance at -50 cm H₂O Distal Hydrostatic Pressure


The Strata II valves are fabricated with a precision molded firm polypropylene and acetal plastic base invested in a smooth, flexible silicone elastomer housing with a contoured bottom (Figs. 3 and 4). The base contributes to structural integrity and pressure/flow reliability, and is designed to prevent distortion of the valve mechanism, siphon control ports and diaphragms. The dome is marked with a radiopaque tantalum-impregnated direction arrow showing flow direction. The currently set Performance Level is discernible through use of the Strata II Adjustment Tools, the StrataVarius Adjustment System or by X-ray (for corresponding settings, refer to Fig. 5).



The Strata II valves incorporate connectors designed to decrease the possibility of catheter disconnection. One encircling ligature is required to secure each catheter to the valve. Radiopaque markers at the base of each connector (Figs. 3 and 4) allow X-ray visualization of the relative positions of catheters and valve in vivo when radiopaque catheters are used.

The dome is designed to allow injection or CSF sampling via a 25-gauge or smaller noncoring needle. Complete valve penetration by the needle is prevented by the firm plastic base. The Strata II valve can be flushed in either the distal or proximal direction by percutaneous finger pressure. Occluders are located proximal and distal to the central reservoir on the Strata II valve regular and distal only on the small model to facilitate selective flushing.

3 Strata II Valve Regular

Regular

4 Strata II Valve Small**5****Patient Data Survey System**

Each Strata II valve is packaged with a patient wallet card and a set of adhesive-backed labels. One label **for each shunt component** can be put on the chart in the operating room, on the patient's clinic chart, and on the patient card. The card can be carried by the patient at all times, providing a current record of implanted devices. Be sure to record the performance level setting on the patient card during the initial setting, and for subsequent revisions to the performance level setting.

Indications

The Strata II valve is a shunt component designed to provide continued cerebrospinal fluid (CSF) flow from the ventricles of the brain into the right atrium of the heart or the peritoneal cavity. The Strata II valve allows the physician to non-invasively adjust the pressure/flow performance level pre- and post-implantation without the need for radiographic confirmation in order to address changing patient needs.

Contraindications

Shunting of CSF into the peritoneal cavity or other areas of the body should not be carried out if there is infection in any areas in which the various components of the shunt system will be implanted. These include infections of the scalp and other skin area through which the shunt system will traverse, the meninges and cerebral ventricles, peritoneum and intraperitoneal and retroperitoneal organs, pleura and blood stream. CSF shunting is contraindicated if there is infection present in any area of the body. Additionally, shunting into the atrium of patients with congenital heart disease or other serious cardiopulmonary abnormalities is contraindicated.

Warnings and Precautions

The valve pressure level setting should always be verified following patient exposure to high magnetic fields.

Devices known to contain magnets should be kept away from the immediate valve implant location, as they may have an effect on the performance level setting of the Strata-type valve. All magnets have an exponentially decreasing effect on the valve the greater the distance from the valve.

Medtronic recommend that all products with magnets used in daily lives to be kept a minimum of 2" (5cm) away from the site where the valve is implanted. Common environmental levels of electromagnetic (radio frequency) radiation generated by Bluetooth, computers, security scanners, metal detectors, microwave ovens, mobile telephones, high voltage lines, and transformers should not affect the performance level settings.

Walk-through security detectors in stores and airports, as well as hand-held security wands, do not have enough magnetic strength to change the setting of the adjustable valve. The full body imaging systems used at airports are x-ray based and will not interfere with the valve setting. Bluetooth is a registered trademark of Bluetooth SIG, Inc.

Valve function and performance level setting should be checked in the event that the valve is subjected to significant mechanical shock or trauma.

Use the Medtronic Neurosurgery Adjustment Kit **REF 45805** or the StrataVarius Adjustment System **REF 45806** to change Performance Level in the Strata II valve.

The Adjustment Tool contains strong magnets. Care should be taken when using the tool near magnetically sensitive medical implants (e.g. pacemakers and vagal nerve stimulators), electronic equipment, data storage devices such as computer diskettes or credit cards.

The Strata II Adjustment Kit and the StrataVarius Adjustment System should **NOT** be sterilized.

Ferromagnetic substances may impede the ability of the adjustment tools to change and confirm the Performance Level setting.

Refer to the Instructions for Use which accompany the Strata II Adjustment Kit or the StrataVarius Adjustment System for instructions, warnings, precautions and complications.

The appropriate product and size must be chosen for the specific patient's needs, based on diagnostic tests and physician experience. Product labeling specifies applicable product performance levels or ranges.

Lint, fingerprints, talc, other surface contaminants, or residues from latex gloves can cause foreign body or allergic reactions.

Improper use of instruments in the handling or implantation of shunt products may result in the cutting, slitting or crushing of components. Such damage may lead to loss of shunt integrity, and necessitate premature surgical revision of the shunt system.

Care must be taken to ensure that particulate contaminants are not introduced into shunt components during preimplantation testing or handling. Introduction of contaminants could result in improper performance of the shunt system. Particulate matter that enters the shunt system may result in shunt occlusion, or may also hold pressure/flow controlling mechanisms open, resulting in overdrainage.

In securing catheters to connectors, the encircling ligatures should be securely, but not too tightly, fastened, lest they eventually cut through the silicone tubing.

Care must be taken in the routing of catheters to prevent kinking and needless abrasion along their course. Abrasion can result in premature catheter failure (fracture). The rim of the twist drill or burr hole may be trimmed to provide a beveled notch where the ventricular catheter emerges and is curved to lie adjacent to the skull.

"Small" size catheters have thinner walls and lower overall strength as compared with "Standard" size catheters. These characteristics result in a comparatively greater potential failure (fracture) rate and, therefore, shorter life expectancy for "Small" size catheters. Physicians who implant "Small" size catheters for cosmetic reasons must acknowledge the potentially higher rate of catheter revision and weight this against the cosmetic benefit. Patients with hydrocephalus shunt systems must be kept under close observation in the postoperative period for signs and symptoms that suggest shunt malfunction. The clinical findings may indicate shunt malfunction. The clinical findings may indicate shunt obstruction or overdrainage of CSF. Shunt obstruction may occur in any of the components of the shunt system. The system may become occluded internally due to tissue fragments, blood clots, tumor cell aggregates, bacterial colonization or other debris. Catheters which contact internal body structures can become kinked or blocked at their tips (e.g., investment of a ventricular catheter tip into the choroid plexus or of the distal catheter tip into the greater omentum or loops of the bowel). Finally, shunt obstruction may occur due to growth of an infant or child, or physical activities which result in disconnection of the shunt components or withdrawal of a distal catheter from its intended drainage site.

Shunt obstruction may occur in any of the components of the shunt system. The ventricular catheter may become occluded by particulate matter such as blood clots or brain fragments, by investment of the catheter tip in choroid plexus, by embedding of the catheter in brain tissue, or by coaptation of the ventricular walls in the presence of overdrainage ("slit ventricle").

Disconnected shunt components may further migrate.

Shunt systems may fail due to mechanical malfunction (example: improper connection of system), leading to under- or overdrainage.

Malfunction or obstruction of the shunt system may lead to signs and symptoms of increased intracranial pressure if the hydrocephalus is not compensated. In the infant, the common symptoms are increased tension of the anterior fontanelle, congestion of scalp veins, listlessness, drowsiness and irritability, vomiting and nuchal rigidity. In older children and adults, the common symptoms are headaches, vomiting, blurring of vision, nuchal rigidity, deterioration of consciousness, and various abnormal neurological findings.

Overdrainage of CSF may predispose development of a subdural hematoma or hydroma or collapse of the lateral ventricular walls leading to obstruction of the ventricular catheter.

If the ventricular catheter becomes bound to the choroid plexus or adjacent brain tissue by fibrous tissue adhesions, it is suggested that it should not be forcibly removed. It is suggested that **gentle** rotation of the catheter may help to free it. It is advised that the catheter be left in place rather than risk intraventricular hemorrhage which may be caused by forcible removal.

Subcutaneous catheter passers can break at welds or component assembly points, or due to extreme deformation of the malleable shaft. Sudden breakage can lead to trauma of tissues or organs, and damage to the shunt system. Instruments must be inspected prior to use to ensure continued integrity and functionality. Disposable instruments must never be reused, or injury to the patient and physician is possible.

Clotting around the atrial portion of a catheter may lead to embolization of the pulmonary arterial tree with resulting cor pulmonale and pulmonary hypertension.

An implanted valve may interfere with a diagnostic MRI if the area in question is near the location of the valve.

Risks associated with radiographic imaging include exposure to low levels of radiation, radiation-induced injuries to the skin and underlying tissues, or the possibility of developing a radiation-induced cancer. Care should be taken to use the least amount of radiation exposure needed to produce the image.

The performance characteristics of this device may be altered if components or devices of other manufacturers are used in conjunction with this device. Testing shows that applying a water column of up to 100cm positive pressure or 50cm negative pressure to the shunt system does not impair performance.

The performance characteristics are verified for all Medtronic shunt components, which are designed to be interchangeably connected with each other. The use of components or devices of other manufacturers in conjunction with this device has not been verified.

MRI Information

Valves

The Strata II valve is considered Magnetic Resonance Conditional in accordance with ASTM F2503.

MRI systems of up to **3.0 Tesla** may be used any time after implantation and will not damage the Strata II valve mechanism, but can change the performance level setting. The performance level setting should always be checked before and after MRI exposure.

The results of the tests performed to assess magnetic field interactions, artifacts, and heating, indicated the presence of the valves evaluated should present no substantial risk to a patient undergoing an MRI procedure using the following conditions:

- Static magnetic field of 3.0 Tesla or less
- Spatial Gradient of 720 G/cm or less
- Radio Frequency (RF) Fields with an average Specific Absorption Rate (SAR) of 3 W/kg for 15 minutes.

Using the GE 3.0T Excite™ HD Magnetic Resonance Imaging System, the valve experienced a maximum temperature change of 0.4°C over a 15-minute exposure period.

The table provides maximum signal voids (artifact sizes) for standard imaging pulse sequences at **3.0 Tesla** per ASTM F2119.

Valve	Pulse Sequence	Plane Imaging	Max. Signal Void (Artifact), cm ²
Strata II	T1-SE	Parallel	35.16
	T1-SE	Perpendicular	33.03
	GRE	Parallel	73.94
	GRE	Perpendicular	66.55

Biological debris inside the valve may impact adjustability, and may lead to adjustment mechanism damage if exposed to 3.0 Tesla MRI. If difficulty is experienced adjusting or reading the valve setting, radiographic setting confirmation should be considered. The reading from the Strata II Indicator tool or the StrataVarius system may be reversed (180 degrees opposite) from the radiographic image. In this situation, radiographic imaging should be used to determine the setting of the valve.

Strata II Adjustment Kit and StrataVarius Adjustment System

Do **NOT** take the Adjustment Tools or the StrataVarius system into an MRI facility as the magnets could potentially be a safety hazard to the patient and/or user.

Proximity to MRI suite may impede the verification mechanisms due to the field strength of an MRI magnet. Move out of the vicinity prior to attempting to verify a valve setting.

Instructions for Use

Setting the Performance Level Before Implantation

The physician must determine the proper initial performance level setting for each patient and adjust the valve accordingly prior to implantation.

Strata II valve PL settings are adjusted with the Strata II Adjustment Kit [REF] 45805 (Fig. 6) or the StrataVarius Adjustment System [REF] 45806.

Refer to the Instructions for Use, which accompany the Adjustment Tools or the StrataVarius Adjustment System for information regarding Performance Level adjustments prior to implantation.

Note: The Strata II Adjustment Tool is labeled [REF] 45805 (Fig. 6).

WARNING: STRATA II VALVES MAY NOT BE RESTERILIZED. IT IS ESSENTIAL TO SET THE PRESSURE LEVEL BEFORE IMPLANTATION WITH THE VALVE SEALED IN ITS STERILE PACKAGE.

Patency Check

Place the inlet connector of the valve into filtered, sterile isotonic saline.

Depress and release the valve dome repeatedly until fluid flows from the outlet connector (Fig. 7). If fluid flows from the outlet connector each time the dome is depressed, the valve is patent.

CAUTION: DO NOT ATTEMPT TO ASPIRATE FLUID THROUGH THE STRATA II VALVE FROM THE DISTAL (OUTLET) END. DAMAGE TO THE VALVE MAY RESULT.


CAUTION: PARTICULATE MATTER IN SOLUTIONS USED TO TEST VALVES MAY RESULT IN IMPROPER PRODUCT PERFORMANCE.

Optional Preimplantation Performance Level Test

Medtronic Neurosurgery does not recommend functional testing of the valve prior to implantation, due to additional handling and potential contamination. If testing is considered necessary, instructions are available upon request.

Surgical Technique

A variety of surgical techniques may be used in placing the Strata II Valve. The valve is implanted with the flat surface adjacent to the pericranium. Site of placement is at the discretion of the surgeon. To reduce the possibility of post-operative valve migration (e.g., as the result of magnetic influence due to MRI), suture the valve to adjacent tissue by passing a suture through the polyester-fabric-reinforced flanges.

6

7

Placement of a Strata II Valve will have an impact on overall shunt performance (Fig. 8). The foramen of Monro and the tip of the ventricular catheter are established zero level reference points for intracranial pressure monitoring. Placement of the valve above these reference points will result in an overall increase in the resistance to flow of the shunt system, with potential underdrainage. Placement of the valve below these reference points will result in an overall decrease in the resistance to flow of the shunt system, with potential overdrainage. These factors must be taken into account for proper patient therapy.

CAUTION: VALVE PLACEMENT CAN HAVE SIGNIFICANT EFFECTS ON THE PERFORMANCE OF THE STRATA II VALVE.

The resistance properties of Strata II Valves as indicated in the Performance Characteristic charts in this product labeling are specified at the zero reference point.

It is suggested that the Strata II valve be placed in a surgically created loose subgaleal pocket, avoiding compression by the overlying scalp, and not under the scalp incision. The Strata II valve should not be placed under the skin of the neck, chest, or abdomen.

Connect valve to catheters by inserting integral connectors into catheters. The connectors should be completely covered by catheter tubing. Secure catheters to connectors with encircling ligatures.

All Medtronic CSF ventricular catheters are compatible w/ the proximal end of Strata II valves. For the distal valve end, all Medtronic CSF cardiac/peritoneal and CSF peritoneal catheters are compatible with Strata II valves. The Small Lumen Peritoneal Catheter is not recommended for use with any Medtronic valve that includes an integrated Delta Chamber mechanism as it affects the overall expected shunt performance. The use of components or devices of other manufacturers in conjunction with Medtronic valves has not been verified.

Injection into the Valve

The Strata II valve is designed to allow injection or CSF sampling through the dome by use of a 25-gauge or smaller noncoring needle (Fig. 9).

The needle should be inserted at an angle of approximately 30° to 45° from the scalp or base of the valve. If the valve will be punctured several times, it is recommended that the needle be inserted at various locations to avoid multiple punctures at a single point. The catheter tubing, occluders, and siphon control diaphragms should **not** be used as injection sites.

CAUTION: LOW TEAR STRENGTH IS A CHARACTERISTIC OF MOST SILICONE ELASTOMER MATERIALS. CARE MUST BE TAKEN ON INSERTION AND REMOVAL OF THE NEEDLE. IMPROPER INSERTION AND/OR REMOVAL OF THE NEEDLE MAY LEAD TO DEVICE LEAKAGE.

CAUTION: THE SIPHON CONTROL PORTION OF THIS DEVICE IS NOT DESIGNED TO ALLOW NEEDLE PENETRATION. PUNCTURE OF THE DIAPHRAGMS MAY AFFECT THE PERFORMANCE CHARACTERISTICS OF THE VALVE AND COMPROMISE THE SHUNT SYSTEM.

Valve Flushing

The Strata II valve can be flushed in both the proximal and distal directions by percutaneous depression of the valve dome (Fig. 10). To selectively flush the Strata II valve regular, depress and occlude either proximal or distal occluder sections of the valve by percutaneous finger pressure, then depress the valve dome. For the small valve, distal flushing can be achieved by depressing the ventricular catheter tubing for occlusion.

The selective flushing will cause fluid to flow in the direction opposite the occluded side of the valve. If there is noticeable resistance to compression, the catheter being flushed may be occluded.




CAUTION: SHUNT OBSTRUCTION MAY OCCUR IN ANY COMPONENT OF A SHUNT SYSTEM AND SHOULD BE DIAGNOSED BY CLINICAL FINDINGS AND DIAGNOSTIC TESTING. VALVE FLUSHING CHARACTERISTICS MAY NOT BE ADEQUATE TO DIAGNOSE OCCLUSION OF CATHETERS. SEE WARNING SECTION.

Setting the Performance Level After Implantation

Strata II valve PL settings are adjusted with the Strata II Adjustment Kit [REF] 45805 or the StrataVarius Adjustment System [REF] 45806. Refer to the **Instructions for Use** which accompanies the Adjustment tools or the StrataVarius Adjustment System for information regarding Performance Level adjustments after implantation.

Note: The Strata II Adjustment Tool is labeled [REF] 45805 (Fig. 6).

CAUTION: EXCESSIVE SWELLING OR BANDAGES MAY MAKE IT DIFFICULT TO DETERMINE A VALVE SETTING. WAIT UNTIL SWELLING IS REDUCED, OR CONFIRM WITH X-RAY. IN ADDITION, SCALP THICKNESS CAN MAKE IT DIFFICULT TO DETERMINE VALVE SETTING. IF THIS IS THE CASE, CONFIRM VALVE SETTING WITH X-RAY.

How Supplied

Medtronic Neurosurgery Strata II valves are packaged **STERILE** and **NON-PYROGENIC** and are intended for **single** (one-time) **use only**.

The enclosed product is designed for single patient use only. Do not re-use, re-process, or re-sterilize this product. Re-use, re-processing or re-sterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death. Do not use if package has been previously opened or damaged. Medtronic Neurosurgery is not responsible for the performance of any product that has been resterilized.

A complete shunt requires a Ventricular Catheter, a Strata II valve, and a Distal Catheter.

Storage Condition

Store devices in dry, clean conditions at normal indoor ambient room temperature.

Patient Education

It is the responsibility of the physician to educate the patient and/or their representative(s) regarding CSF shunting. This should include a discussion of warnings and precautions highlighted in these instructions, a description of the complications associated with implantable shunt systems, and an explanation of potential alternative products and treatments.

Complications

Complications associated with ventriculoperitoneal CSF shunting systems may be similar to those experienced in any surgical procedure carried out under local and/or general anesthesia. These include reactions to drugs and anesthetic agents, electrolyte imbalance and excessive blood loss, particularly in infants. A patient may rarely exhibit a reaction due to sensitivity to the implant.

In CSF shunting procedures, the most common complications are due to obstruction of the system as described under "Warnings." Obstruction may occur in any component of the system due to plugging by brain fragments, blood clots, and/or tumor cell aggregates at some point along its course. Obstruction may also occur because of separation of the system components or kinking and/or coiling of the catheter. This may predispose migration of the ventricular catheter into the lateral ventricle and the distal catheter into the peritoneum, or other structure in which the catheter is implanted. As noted previously, growth of the infant or child may cause the distal catheter to be withdrawn from the atrium into the internal jugular vein or from the peritoneum into tissue planes where the fluid cannot be absorbed.

There are other potentially serious complications. Local and systemic infections are not uncommon with shunting procedures. Usually, they are due to organisms inhabiting the skin, particularly *Staphylococcus epidermidis*. Other pathogens circulating in the blood stream may colonize the shunt and, in the majority of patients, require its removal.

In 1973, Robertson et al. summarized the incidence of infection in ventriculoperitoneal shunts reported up to that time. Infection in ventriculoperitoneal shunting occurred in 5 to 10% of the patients in most of the reports.

In 1993, Kettle et al. reported significant reductions in infection (less than 4%) with the use of antibiotics, short duration of surgery (surgical experience) and control of the operating room environment (e.g., designated operating room, limited personnel and traffic, covered skin surfaces). The article states that results can also be obtained without the use of antibiotics, but with rigorous perioperative control of the environment.

Using prophylactic antibiotics in shunted patients is somewhat controversial as their use may predispose infection by more resistant organisms. Therefore, the decision to use antibiotics prophylactically rests with the attending physician and/or surgeon.

Shunting into the peritoneum may fail because of investments of the catheter in loops of bowel or in the greater omentum. Perforation of the bowel by the peritoneal catheter with subsequent development of peritonitis has been described.

CSF overdrainage may result in excessive reduction of CSF pressure and predispose the development of a subdural hematoma or hygroma, and excessive reduction of ventricular size leading to obstruction because of impingement of the ventricular walls on the inlet holes in the catheter. In the infant, this excessive pressure reduction will cause marked depression of the anterior fontanelle, overriding of cranial bones and may convert communicating into obstructive hydrocephalus.

The incidence of epilepsy after ventricular shunting procedures has been reported. This study also indicated that the incidence of seizures increased with multiple catheter revisions.

Intra-abdominal complications associated with peritoneal shunting include perforation of the small or large bowel with resultant peritonitis, perforation of other viscera, ureteral obstruction, migrations of the catheter into unintended anatomical location, and development of ascites and pseudocysts.

Returned Goods Policy

Products must be returned in unopened packages, with manufacturer's seals intact, to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling. Determination of a product defect or mislabeling will be made by Medtronic Neurosurgery, which determination will be final. Products will not be accepted for replacement or credit if they have been in the possession of the customer for more than 90 days.

Warranty

A. Standard Limited Warranty. Medtronic Neurosurgery warrants to the original end user purchaser ("Purchaser") that Medtronic Neurosurgery's single-use implantable product ("Product") purchased by Purchaser, at the time of delivery to Purchaser, shall be substantially free from defects in material and workmanship. Medtronic Neurosurgery makes no warranty (express, implied, or statutory) for Products that are modified (except as expressly contemplated herein) or subjected to unusual physical stress, misuse, improper operation, neglect, improper testing, use in combination with other products or components other than those for which the Products were designed, or use in any manner or medical procedure for which the Products are not indicated.

B. Remedy. Purchaser's exclusive remedy and Medtronic Neurosurgery's sole liability for breach of the foregoing warranty shall be, at Medtronic Neurosurgery's sole option and election, to replace the Product or credit Purchaser for the net amount actually paid for any such Product; provided that (i) Medtronic Neurosurgery is notified in writing within ninety (90) days after Purchaser's receipt of the Product that such Product failed to conform, including a detailed explanation in English of any alleged nonconformity; (ii) such Product is returned to Medtronic Neurosurgery within ninety (90) days after Purchaser's receipt of the Product F.O.B. 125 Cremona Drive, Goleta, California 93117, U.S.A. or as otherwise designated by Medtronic Neurosurgery; and (iii) Medtronic Neurosurgery is reasonably satisfied that the claimed nonconformities actually exist. Except as expressly provided in this paragraph, Purchaser shall not have the right to return Products to Medtronic Neurosurgery without Medtronic Neurosurgery's prior written consent.

C. Exclusion of Other Warranties. EXCEPT FOR THE LIMITED WARRANTY PROVIDED IN (A) ABOVE, MEDTRONIC NEUROSURGERY GRANTS NO OTHER WARRANTIES OR CONDITIONS, EXPRESS OR IMPLIED, AND MANUFACTURER SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTIES AND CONDITIONS OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC NEUROSURGERY NEITHER ASSUMES NOR AUTHORIZES ANY OTHER PERSON TO ASSUME ANY OTHER LIABILITIES ARISING OUT OF OR IN CONNECTION WITH THE SALE OR USE OF ANY PRODUCT.

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© 2018 Medtronic, Inc.
Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic, Inc.