Medtronic

CG FUTURE[®] 638R/638B

Annuloplasty Ring/Band - Mitral

Annuloplastikkring/-bånd - mitral

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1	DEVICE DESCRIPTION	5
2	INDICATIONS	5
3	CONTRAINDICATIONS	5
4	WARNINGS	5
5	PRECAUTIONS	6
6	POTENTIAL ADVERSE EVENTS	6
7	INDIVIDUALIZATION OF TREATMENT	7
8	PATIENT COUNSELING INFORMATION	7
9	HOW SUPPLIED	7
	9.1 Packaging	7
	9.2 Storage	7
10	DIRECTIONS FOR USE	8
	10.1 Sizing	8
	10.2 Handling and Preparation Instructions	8
	10.3 Device Implantation	8
	10.4 Accessories	10
	10.5 Sterilization	10
11	REGISTRATION INFORMATION	10
12	DISCLAIMER OF WARRANTIES	11

1 DEVICE DESCRIPTION

The Medtronic CG Future® Annuloplasty Ring Model 638R and the Medtronic CG Future® Annuloplasty Band Model 638B (hereafter referred to as CG Future® Annuloplasty Ring/Band) consist of a MP35N wire overmolded with silicone and covered with polyester fabric. Both the ring and band must be implanted in the mitral position. The ring and band are marked at three points by colored sutures; two markers correspond to the trigones of the mitral valve and one identifies the midpoint of the device. The device size is identified by the inside diameter of the ring or band at its widest point. The MP35N wire allows radiographic visualization of the devices. The ring also includes a silicone marker, impregnated with barium sulfate, for radiographic visualization.

Nonclinical testing has demonstrated that the CG Future® Annuloplasty Ring/Band is MR (magnetic resonance) Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of ≤ 3.0 tesla
- Spatial gradient field of ≤ 3.9 T/m
- Maximum whole-body-averaged specific absorption rate (SAR) of 1.1 W/kg for ≤ 20 minutes of scanning as read from equipment monitor

In nonclinical testing the CG Future® Annuloplasty Ring/Band produced a temperature rise of $\leq 0.2^{\circ}\text{C}$ at a maximum whole-body-averaged specific absorption rate (SAR) of 1.1 W/kg for ≤ 20 minutes of MR scanning in a field strength of 3.0 tesla, GE Signa®¹ LX 3.0 T MR System: multi-coil superconducting, actively shielded magnet housed in a 94 cm bore; gradient amplitude - 40 mT/m and slew rate - 150 T/m/s; VH3_M4 software. The maximum whole body SAR of 1.1 W/kg for ≤ 20 minutes was read from the scanner console.

2 INDICATIONS

The CG Future® Annuloplasty Ring/Band are indicated for the reconstruction and/or remodeling of pathological mitral valves. Valvular insufficiency and/or stenosis may be corrected by appropriate repair and annular remodeling.

3 CONTRAINDICATIONS

- Heavily calcified valves.
- Valvular retraction with severely reduced mobility.
- Active bacterial endocarditis.

4 WARNINGS

- **For single use only.**
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring and band implant and sizing techniques, should use this device.

¹ Signa® is a registered trademark of General Electric Corporation.

- Correct annuloplasty ring and band sizing is an important element of a successful valve repair. Undersizing the ring or band can result in valve stenosis, ring/band dehiscence, and/or ring/band fracture. Oversizing the ring or band can result in valve regurgitation and/or ring/band fracture.
- Care should be taken that sutures are not placed in atrial tissue, as this may result in impairment of the cardiac conduction system.
- It is necessary to secure the ends of the ring and band into fibrous trigone tissue to maintain permanent attachment.
- Suture knots must be securely tied. Loose knots and long suture tails may be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the ring or band, as resultant loose threads may be a source of hemolysis, thrombosis, and/or thromboembolism.
- Intraoperative and/or postoperative echocardiography should be used to evaluate the effectiveness of the valve repair. Minimizing regurgitation and systolic anterior motion (SAM) are important elements of an effective repair.
- Surgeons who use annuloplasty rings and bands should be current on all anticoagulation regimens.
- When postoperative anticoagulant therapy is used, the patient's anticoagulation status should be carefully monitored.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulation therapy.
- The surgeon may desire that patients in atrial fibrillation remain on anticoagulation therapy until sinus rhythm is established.
- Do not alter or deform the ring or band to conform to annular anatomy as this could lead to ring/band fracture and possible mitral regurgitation.
- Do not squeeze the ring or band with sharp instruments as this may damage the surface of the stiffener, which may result in ring/band fracture and possible mitral regurgitation.
- The ring and band are indicated for use in the mitral position only. Use in the tricuspid position may result in ring/band fracture, conduction system damage, tricuspid regurgitation, and/or ring/band dehiscence.

5 PRECAUTIONS

- Do not use cutting edge needles, as they may damage the annuloplasty device potentially leading to ring/band dehiscence, ring/band fracture, and possible mitral regurgitation.

6 POTENTIAL ADVERSE EVENTS

While infrequent, certain complications have been reported when using annuloplasty ring/bands. These include the following:

- Uncorrected or recurrent regurgitation
- Stenosis
- Ring/Band dehiscence
- Hemolysis (even with mild regurgitation)

- Low cardiac output
- Heart block
- Systolic anterior motion (SAM) and left ventricular outflow tract obstruction (LVOTO)
- Damage to coronary arteries
- Endocarditis
- Thrombosis
- Thromboembolism
- Anticoagulant-related hemorrhage
- Ring/Band fracture
- Leaflet perforation

The potential for these complications should be considered when selecting the most beneficial surgical procedure for each patient.

To avoid or minimize occurrence of these adverse events, the annuloplasty repair, including sizing and implantation, should be conducted in accordance with the methods described in these Instructions for Use by surgeons with appropriate training and experience in valve repair.

7 INDIVIDUALIZATION OF TREATMENT

To allow for healing and incorporation of the annuloplasty ring or band by host tissue, regardless of cardiac rhythm, postoperative anticoagulation therapy should be considered for at least six weeks following surgery.

8 PATIENT COUNSELING INFORMATION

Patients with annuloplasty rings or bands who undergo dental or other potentially bacteremic procedures must be considered for prophylactic antibiotic therapy.

9 HOW SUPPLIED

9.1 Packaging

The CG Future® 638R Annuloplasty Ring is available in the following sizes for the mitral position: 24, 26, 28, 30, 32, 34, 36, and 38 mm. The CG Future® 638B Annuloplasty Band is available in the following sizes: 26, 28, 30, 32, 34, 36, and 38 mm. The package contains a single annuloplasty ring or band assembly consisting of the device and holder (Figure 1) packaged in sterile, double-aseptic transfer pouches. The packaging system is designed to ease placement of the device into the sterile field. The ring and band assemblies are sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are NONSTERILE and must not be placed in the sterile field.

9.2 Storage

Store the product in its original packaging, including the outer shelf carton, in a clean, cool, and dry area to protect the product and minimize the potential for contamination.

The sterility and nonpyrogenicity of the CG Future® Annuloplasty Ring/Band are validated to remain unaffected until the Use By date identified on the shelf carton, provided the pouches are not opened or damaged.

10 DIRECTIONS FOR USE

10.1 Sizing

The CG Future® 7638 Ring/Band Sizer Set is used to select the proper CG Future® Annuloplasty Ring/Band size for repair of the mitral valve. Sizers are reusable; however, they must be cleaned and sterilized by autoclave (steam) prior to each use.

Proper ring/band size selection is an important part of valvular annuloplasty to help restore proper function. Use the CG Future® 7638 Annuloplasty Ring/Band Sizer Set and Medtronic 7615 Annuloplasty Handle for size selection.

The handle and sizers are provided NONSTERILE; they must be cleaned and sterilized prior to use. Align and insert the handle into the snap cavity of the sizer (Figure 2). The thin section of the handle is malleable, allowing the surgeon to align the sizers with the valve annulus.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the CG Future® Annuloplasty Ring/Band.

Mitral Valve

To determine the proper ring or band size, both the distance between the annular trigones and the area of the anterior leaflet must be measured. First, lower the sizer onto the valve annulus and align the sizer notches with the annular trigones. Second, gently extend the anterior (A) leaflet and cover its surface with the selected sizer. The sizer that has a notch spacing most nearly matching the intertrigonal distance (Figure 3) and a surface area most nearly matching that of the anterior leaflet corresponds to the size of the ring or band that should be selected.

10.2 Handling and Preparation Instructions

- Open the box and remove the product literature and Patient Registration Form.
- Remove the double-aseptic transfer pouch containing the device mounted to the holder.
- Inspect the pouches ensuring they have not been opened or damaged. The ring and band assemblies are sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the device.
- If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
- Open the outer transfer pouch, and while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
- The inner pouch should be opened only in the sterile field.

10.3 Device Implantation

The CG Future® Annuloplasty Ring/Band assembly may be used with or without the handle (Medtronic 7615 Annuloplasty Handle, provided separately).

Caution: The CG Future® Annuloplasty Ring/Band are designed for implantation only with interrupted-suture techniques.

Warning: The handle must be cleaned and sterilized prior to use. To use the handle, align and engage the sterile handle into the snap fit cavity of the holder (Figure 4). The thin section of the handle is malleable, allowing the surgeon to align the device with the annulus.

Note: The handle must be inserted into the side of the holder with the snap fit cavity and the printed text.

Remove the serial number identification tag (Figure 5) and record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

Warning: The serial tag must be removed from the ring/band for proper function. Do not cut or tear the ring/band fabric during removal of the serial tag.

Mitral Annulus Suture Placement

Place sutures in each trigone, approximately 4 mm in width. Place additional interrupted sutures, approximately 4 mm in width, in the posterior portion of the mitral annulus following the line of cusp insertion (Figure 6). A total of approximately 10 to 14 sutures should be placed for the ring/band.

Warning: Avoid placing sutures in the circumflex coronary artery.

Note: Pledgets may be used for trigone sutures and mid-posterior annulus suture to reduce the possibility of ring/band dehiscence.

Band Suture Placement

Orient the band assembly on the mitral (M) annulus on the holder (Figure 7).

Pass sutures through the band, approximately 2 mm in width, entering at the bottom of the band and exiting the periphery of the band (Figure 8). The sutures placed in the trigonal tissue should be passed through the eyelets of the band stiffener (Figure 9). Pass only one suture through the trigonal marker on the device.

Ring Suture Placement

Orient the ring assembly on the mitral (M) annulus on the holder (Figure 7).

Pass sutures through the ring, approximately 2 mm in width, entering at the bottom of the ring and exiting the periphery of the ring (Figure 8). Place sutures completely around the device, including the anterior mitral annulus region. The sutures placed in the trigonal tissue should be passed through the eyelets of the ring stiffener (Figure 9). Pass only one suture through the trigonal marker on the device.

Caution: The ring/band must be sutured through the eyelets into the trigonal tissue to ensure secure attachment.

Holder/Handle Removal

Use the holder/handle to push the ring or band down onto the valve annulus while pulling back on the sutures. Using a scalpel, cut the retention sutures on the holder in the areas indicated (Figure 10). Remove the disposable holder and handle from the valve annulus (Figure 11). Dispose of the holder.

Warning: Do not cut the ring or band fabric while cutting the holder retention sutures.

Knots, Testing

Tie all knots around the device securely (Figure 12), trim all excess sutures, and test valvular competency after removing the holder.

Warning: The CG Future® Annuloplasty Ring/Band holder must be removed from the device at the end of the procedure for proper function. Under no circumstances should the holder remain attached to the ring or band.

Warning: Suture knots must be securely tied. Loose knots and long suture tails may be a source of hemolysis, thrombosis, or thromboembolism.

10.4 Accessories

Use the CG Future® 7638 Annuloplasty Ring/Band Sizer Set to determine the appropriate device size. Refer to the CG Future® Annuloplasty Ring/Band Sizer Set Instructions for Use for detailed information on accessory use and sterilization.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the CG Future® Annuloplasty Ring/Band.

Use only the Medtronic 7615 Annuloplasty Handle to interface with the holder and sizers. Refer to the Medtronic Annuloplasty Handle Instructions for Use for detailed information on use and sterilization.

10.5 Sterilization

10.5.1 Ring/Band (Model 638R/638B)

The ring and band are provided sterile (ethylene oxide) on the holder and must not be resterilized. Rings and bands that have been damaged or contaminated from patient contact should not be used.

11 REGISTRATION INFORMATION

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU.

A patient registration form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the device. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

12 DISCLAIMER OF WARRANTIES

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

DISCLAIMER OF WARRANTY

ALTHOUGH THE CG FUTURE® 638R/638B ANNULOPLASTY RING/BAND, THE CG FUTURE® 7638 RING/BAND SIZER SET, AND THE MEDTRONIC 7615 ANNULOPLASTY HANDLE, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY

ALTHOUGH THE CG FUTURE® 638R/638B ANNULOPLASTY RING/BAND, THE CG FUTURE® 7638 RING/BAND SIZER SET, AND THE MEDTRONIC 7615 ANNULOPLASTY HANDLE, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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
1. Antunes MJ, Colsen PR, Kinsley RH. Mitral valvuloplasty: A learning curve. *Circulation*. 1983;68(suppl II):70-75.
2. Carpentier A. Cardiac valve surgery—The French correction. *J Thorac Cardiovasc Surg*. 1983;86:323-337.
3. Carpentier A. Mitral Valve Reconstructive Surgery. In: Robb C, Smith R, eds. *Operative Surgery*. 3rd ed. London, England: Butterworth & Co Ltd; 1977:169-177.
4. Carpentier A. Plastic and reconstructive mitral surgery. In Kalamanson D, ed. *The Mitral Valve. A Pluridisciplinary Approach*. Acton, MA: Publishing Sciences Group; 1976:527-540.
5. Carpentier A, Chauvaud S, Fabiani JN, et al. Reconstructive surgery of mitral valve incompetency: A ten-year appraisal. *J Thorac Cardiovasc Surg*. 1980;79(3):338-348.
6. Deloche A, Jebara VA, Relland JY, et al. Valve repair with Carpentier techniques—The second decade. *J Thorac Cardiovasc Surg*. 1990;99(6):990-1002.
7. Galloway AC, Colvin SB, Spencer FC. Mitral Valve Reconstruction. In: Grillo HC, Austen WG, Wilkins EW, Mathisen DJ, Vlahakes GJ, eds. *Current Therapy in Cardiothoracic Surgery*. Philadelphia, PA: BC Decker, Inc; 1989:352-357.
8. Galloway AC, Anderson RV, Grossi EA, Spencer FC, Colvin, SB. Acquired Heart Disease. In: Schwartz SI, Shires T, Spencer FC, Daly JM, Fischer JE, Galloway AC, eds. *Principles of Surgery*. 7th ed. New York, NY: McGraw-Hill; 1999:845-907.
9. Galloway AC, Colvin SB, Baumann FG, et al. Long-term results of mitral valve reconstruction with Carpentier techniques in 148 patients with mitral insufficiency. *Circulation*. 1988;78(suppl I):97-105.
10. Galloway AC, Colvin SB, Baumann FG, Harty S, Spencer FC. Current concepts of mitral valve reconstruction for mitral insufficiency. *Circulation*. 1988;78(5):1087-1098.
11. Galloway AC, Colvin SB, Baumann FG, et al. A Comparison of mitral valve reconstruction with mitral valve replacement: Intermediate-term results. *Ann Thorac Surg*. 1989;47(5):655-662.
12. Galloway AC, Grossi EA, Spencer FC, Colvin SB. Operative therapy for mitral insufficiency from coronary artery disease. *Semin Thorac Cardiovasc Surg*. 1995;7(4):227-232.
13. Grossi EA, Galloway AC, Parish MA, et al. Experience with twenty-eight cases of systolic anterior motion after mitral valve reconstruction by the Carpentier technique. *J Thorac Cardiovasc Surg*. 1992;103(3):466-470.
14. Grossi EA, Galloway AC, LeBoutillier M, et al. Mitral valve repair in the elderly. *J Am Coll Cardiol*. 1994;14(2):465A.
15. Grossi EA, Galloway AC, Steinberg BM, et al. Severe calcification does not affect long-term outcome of mitral valve repair. *Ann Thorac Surg*. 1994;58(3):685-688.

16. Grossi EA, Steinberg BM, LeBoutillier M III, et al. Decreasing incidence of systolic anterior motion after mitral valve reconstruction. *Circulation*. 1994;90(suppl II):195-197.
17. Grossi EA, Galloway AC, LeBoutillier M III, et al. Anterior leaflet procedures during mitral repair do not adversely influence long-term outcome. *J Am Coll Cardiol*. 1995;25(1):134-136.
18. Grossi EA, Galloway AC, LeBoutillier M, et al. Mitral valve repair in the elderly. *Cardiol Elderly*. 1995;3:269-272.
19. Grossi EA, Galloway AC, Miller JS, et al. Valve repair versus replacement for mitral insufficiency: When is a mechanical valve still indicated? *J Thorac Cardiovasc Surg*. 1998;115(2):389-394.
20. Lessana A, Carbone C, Romano M, et al. Mitral valve repair: Results and the decision-making process in reconstruction. *J Thorac Cardiovasc Surg*. 1990;99(4):622-630.
21. Spencer FC, Galloway AC. Recent developments and evolving techniques of mitral valve reconstruction. *Ann Thorac Surg*. 1998;65:307-313.
22. Spencer FC, Galloway AC, Colvin SB. Acquired Diseases of the Mitral Valve. In: Sabiston DC, ed. *Surgery of the Chest*. 6th ed. Philadelphia, PA: WB Saunders Co; 1995:1673-1701.9.
23. Spencer FC, Colvin SB, Culliford AT, Isom OW. Experiences with the Carpentier techniques of mitral valve reconstruction in 103 patients (1980-1985). *J Thorac Cardiovasc Surg*. 1985;90(3):341-350.
24. Virmani R, Chun PK, Parker J, McAllister HA Jr. Suture obliteration of the circumflex coronary artery in three patients undergoing mitral valve operation. *J Thorac Cardiovasc Surg*. 1982;84(5):773-778.

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All Rights Reserved
Medtronic of Canada Ltd.
Medtronic, Inc.
© Medtronic, Inc. 2006, 2007